# Quick PROTOCOL

BIOTEC NSTDA

## "ENZhance Complete" Kit for industrial enzyme detection and screening

### "ENZhance Complete" Components

- 10x ENZhance Reagent (pH 4, 7, 10, filter-sterilized): 40 ml each
- Substrates for cellulase, xylanase, amylase, and protease
- Agar (5 g) (Supplied by users)

- Note:**
- For 10 assays of each enzyme pH
  - Keep ENZhance Complete kit at room temperature and dried conditions

### Preparation of 1x ENZhance overlay agar

(100 ml for 10 plate assays)

- Dissolve 0.7 g agar in 90 ml distilled water (optional: sterilized by autoclaving at 121 °C, 15 min with/ without substrate(s))
- Add 10 ml of 10x ENZhance Reagent and the target enzyme substrate(s)\* when the agar mixture is cooled down to ≈50 °C.

### Protocol for enzyme detection/ screening

1. Cultivate bacterial samples (pure isolates/ mixed culture) on appropriate agar medium under optimal growth conditions (e.g. at 37 °C for 18-48 hrs).
2. Overlay the bacterial culture with 1x ENZhance overlay agar containing the substrate(s) for the target enzyme(s) (≈10 ml).
3. Incubate the plate under enzyme working conditions. Enzyme activity can be detected after 15 min or up to 24 h depending on target enzyme activity level.
4. Detect the target enzyme activity by observing color zones or clear zones around bacterial colonies\*\*.

- Note:**
- Substrate concentration for enzyme detection
    - 0.05% chromogenic substrate for cellulase, xylanase or amylase
    - 2% substrate for protease

\*\* Color zones indicate positive cellulase, xylanase or amylase activity. Clear zones indicate positive protease activity.

## "ENZhance Complete" Kit for industrial enzyme detection and screening

Cultivate bacterial samples (pure isolates/ mixed culture) on appropriate agar medium under optimal growth conditions

Incubation  
(e.g. at 37 °C  
for 18-48 h)

10x ENZhance Reagent  
(pH 4, 7, 10)

Overlay the bacterial culture with 1x ENZhance overlay agar containing the substrate(s) for the target enzyme(s).

Incubate at enzyme working  
temperature 4 - 60 °C  
(15 min - 24 h)

### • Detection of enzyme activity •

Color zones around bacterial colonies indicate positive cellulase, xylanase or amylase activity.

Clear zones around bacterial colonies indicate positive protease activity.

- Note:**
- \*Dissolve 0.7% agar in distilled water (optional: sterilized by autoclaving) and add target enzyme substrate(s) and 10x ENZhance Reagent (1/10 by volume). Overlay when the mixture is cooled down to ≈50 °C. Substrates for other target enzymes of interest are applicable.

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